



Quality by Design (QbD) in Process Analytical Technology (PAT) and Real Time Release Testing (RTRT)

Isranalytica 2012
Tel Aviv, Israel
25 January 2012


Christine M. V. Moore, Ph.D.
Acting Director
ONDQA/CDER/FDA



Outline

- Clarification of QbD, PAT and RTRT Terminology
- Examples of PAT and RTRT
- Regulatory Considerations for RTRT
- Concluding Remarks

New Quality Terminology




What is Quality by Design (QbD)?

- Systematic approach to pharmaceutical development and manufacturing
- Begins with predefined objectives
- Emphasizes product and process understanding and process control
- Based on sound science and quality risk management




From ICH Q8(R2)

Example QbD Approach - ICH Q8(R2)



5

QbD Approach



6



FDA's View of Process Analytical Technologies

- Process Analytical Technology (PAT)
 - a system for designing, analyzing, and controlling manufacturing
 - through timely measurements of critical quality and performance attributes of raw and in-process materials and processes
 - with the goal of ensuring final product quality
- PAT Fundamental Tenets
 - **Quality cannot be tested** into the product; it should be built-in or **should be by design**
- PAT Goals
 - Enhance understanding and control of processes

7



Real Time Release Testing

- Real Time Release Testing (RTRT) is the ability to evaluate and ensure the quality of in-process and/or final product based on process data
 - Typically include a valid combination of measured material attributes and process controls

ICH Q8(R2)


8

Examples of PAT

- During Development
 - In-line laser light scattering analyzer to monitor and understand nucleation during crystallization
 - At-line particle size distribution monitoring
 - NIR to understand & design blending process
- In Manufacturing
 - Table compression weight check and adjustment
 - Endpoint determination of blending
 - Weight check and adjustment of powder filling operation
 - Adjustment of process parameters based on starting material attributes

9


Endpoint Control Example for Blend Uniformity



Uniformity of excipients and blend determined by on-line process monitoring by NIR

10

Example Monitoring and Modeling for Crystallization Control



- Controlled crystallization with seeding
- Processing occurs entirely in thermodynamically favored regime


11

Examples of RTRT Approaches


- On-line or in-line measurements and/or controls, for example
 - Tablet weight after compression
 - Particle size measurement after granulation or milling
 - Moisture measurement during drying
 - Blend uniformity
- Fast at-line measurements, for example
 - NIR for tablet assay
 - Disintegration in lieu of dissolution
- Models as surrogate for traditional release tests, for example
 - Multivariate model as a surrogate for dissolution
- Process signatures
 - *An evolving approach*

12


Example of An Unified Approach for TRT



Conceptual Example of Control Strategy for TRT in Continuous Manufacturing



Multivariate Model for Predicting Dissolution



15

Considerations for Real Time Release Testing Approaches

16

Benefits of RTRT

- Provides for increased assurance of quality
 - More process data collected
- Provides increased manufacturing flexibility and efficiency
 - Shorter cycle time
 - Reduced inventory
 - Reduction in end product testing
 - Reduction in manufacturing cost
- Allows leveraging of enhanced process understanding
 - Corrective actions may be implemented in real time

17

Regulatory Documents Discussing RTRT

- FDA Guidance for Industry: PAT – A Framework for Innovative Pharmaceutical Development, Manufacturing and Quality Assurance”, Sept 2004
<http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM070305.pdf>
- ICH Q8(R2) – Pharmaceutical Development, Aug 2009
http://www.ich.org/fileadmin/Public_Web_Site/ICH_Products/Guidelines/Quality/Q8_R1/Step4/Q8_R2_Guideline.pdf
- ICH Quality Implementation Working Group on Q8, Q9 and Q10, Questions & Answers (R4), Nov 2010
http://www.ich.org/fileadmin/Public_Web_Site/ICH_Products/Guidelines/Quality/Q8_9_10_QAs/Q-1WG_QAs_Step4/Q8_Q9_Q10_Question_and_Answer_R4_step_4_November_2010.pdf
- ICH Quality Implementation Working Group Points to Consider, June 2011
http://www.ich.org/fileadmin/Public_Web_Site/ICH_Products/Guidelines/Quality/Q8_9_10_QAs/PtC/Quality_1WG_PtC_16_June_2011.pdf

18

ICH IWG Q&A's on RTRT

- How is batch release affected by employing RTRT?
- Does RTRT mean elimination of end product testing?
- Is a product specification still necessary in the case of RTRT?
- When using RTRT, is there a need for stability test methods?
- What is the relationship between Control Strategy and RTRT?
- Do traditional sampling approaches apply to RTRT?
- If RTRT results fail or trending toward failure, can end-product testing be used to release the batch?
- What us the relationship between in-process testing and RTRT?
- What is the difference between RTR and RTRT?
- Can surrogate measurement be used for RTRT?

19

ICH IWG Q&A's on RTRT

- How is batch release affected by employing RTRT?
- Does RTRT mean elimination of end product testing?
- Is a product specification still necessary in the case of RTRT?
- When using RTRT, is there a need for stability test methods?
- What is the relationship between Control Strategy and RTRT?
- Do traditional sampling approaches apply to RTRT?
- If RTRT results fail or trending toward failure, can end-product testing be used to release the batch?
- What us the relationship between in-process testing and RTRT?
- What is the difference between RTR and RTRT?
- Can surrogate measurement be used for RTRT?

20

Considerations for the Point of Testing

- Is there a potential for the measured CQA to change downstream from the measurement point? For example,
 - Blend desegregation
 - Loss of weight (e.g., chipping) after weighing step
 - Hydrolytic degradation during aqueous film coating
- Is identity determined at a point that is visually unique?
 - Mitigation of potential human and/or system error
 - Unique identifiers on the intermediate when measured (e.g., embossing, size, shape)
- Risk assessment is valuable to exploring potential failure modes

21

Sampling Considerations

- Probe/sample location representative of entire vessel
- Sample frequency representative of entire batch
- Sample acquisition time
 - Suitable for system dynamics/mixing
- Sample volume/mass
 - Determine amount of sample measured
 - Representative of unit dose
- Sample interface
 - Remains constant over the process (e.g., no fouling)
 - Environmental factors (e.g., temperature, humidity)


22

Considerations for RTRT Specification

- Specification still required in an RTRT approach (CFR §314.50(d) and CFR § 211.165(a))
- Should be representative of actual measurement
 - Can include in-process measurements (e.g., NIR measurements for assay of uncoated tablets)
 - Can include surrogate measurements (e.g., models for dissolution)
 - Methods should be appropriately validated (including models used as surrogate measurements)
- Alternatives can be included for stability monitoring
- Appropriate statistical criteria for large sample sizes

23


Batch Release Decision for RTRT



- In principle, end product testing should not be substituted for failure of an RTRT release method. The failure should be investigated and followed up appropriately.

24

Communicate Early and Often!



25

Interactions with FDA

- Request meeting according to guidance
“Formal Meetings Between the FDA and Sponsors or Applicants”
 - Clearly state as CMC meeting for RTRT
 - Both CMC and GMP questions can be included
- End of Phase II or Pre-supplement submission is a good time to start dialogue
 - Initially, not all details or data are expected to be available
 - Discuss desired or expected approach
 - Ask specific questions
- Additionally, a Pre-Operational Review (or Pre-Operational Visit, POV) can be requested
 - See PAT Guidance (2004) or “ORA Field Management Directive 136”

26



Concluding Remarks

- PAT and RTRT can provide a higher assurance of product quality
 - Real-time analysis and control of process
 - Enhanced process understanding
 - Operational flexibility
 - Framework for continuous manufacturing
 - Support of continual improvement
- ICH has published several documents providing guidance on implementation of RTRT
- FDA supports the implementation of RTRT approaches using a science and risk-based approach
 - Recommend early and frequent discussion with Agency before implementation

27

Thank you!

Questions, comments, concerns:
NewDrugCMC@fda.hhs.gov

28